



Recombinant Bovine Inorganic pyrophosphatase (PPA1)

| | |
|--------------------------|---|
| Product Code | CSB-MP018410BO |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P37980 |
| Product Type | Recombinant Protein |
| Immunogen Species | Bos taurus (Bovine) |
| Purity | >85% (SDS-PAGE) |
| Sequence | SGFSSEERA APFTLEYRVF LKNEKGQYIS PFHDIPIYAD KEVFHMOVVEV PRWSNAKMEI ATKDPLNPIK QDVKKGKRLRY VANLFPYKGY IWNYGAIPT WEDPGHNDKH TGCCGDNDPI DVCEIGSKVC ARGEIIRVKV LGILAMIDEG ETDWKVIAIN VEDPDAANYN DINDVKRLKP GYLEATVDWF RRYKVPDGKP ENEFANAEF KDKNFAIDII ESTHDYWRAL VTKKTDGKGI SCMNNTVSES PFQCDPDAAK AIVDALPPPC ESACTIPTDV DKWFHHQKN |
| Source | Mammalian cell |
| Target Names | PPA1 |
| Protein Names | Recommended name: Inorganic pyrophosphatase EC= 3.6.1.1 Alternative name(s): Pyrophosphate phospho-hydrolase Short name= PPase |
| Expression Region | 2-289 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Target Details | This protein is a member of the inorganic pyrophosphatase (PPase) family. PPases catalyze the hydrolysis of pyrophosphate to inorganic phosphate, which is important for the phosphate metabolism of cells. Studies of a similar protein in bovine suggested a cytoplasmic localization of this enzyme. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |