



# Recombinant Escherichia coli O17:K52:H18 CTP synthase (pyrG)

<b>Product Code</b>	CSB-BP006176EOG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B7N716
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MTTNYIFVTG GVVSSLGKGI AAASLAAILE ARGLNVTIMK LDPYINVDPG          TMSPIQHGEV FVTEDGAETD LDLGHYERFI RTKMSRRNNF TTGRIYSDVL          RKERRGDYLG ATVQVIPHIT NAIKERVLEG GEGHDVVLVE IGGTVGDIES          LPFLEAIRQM AVEIGREHTL FMHLTLVPYM AASGEVKTKP TQHSVKELLS          IGIQPDILIC RSDRAVPANE RAKIALFCNV PEKAVISLKD VDSIYKIPGL          LKSQGLDDYI CKRFSLNCP E ANLSEWEQVI FEEANPVSEV TIGMVGKYIE          LPDAYKSVIE ALKHGGLKNR VSVNIKLIDS QDVETRGVEI LKGLDAILVP          GGFGYRGVEG MITTARFARE NNIPYLGICL GMQVALIDYA RHVANMENAN          STEFVPDCKY PVVALITEWR DENGNEVRS EKSDLGGTMR LGAQQCQLVD          DSLVRQLYNA PTIVERHRHR YEVENMMLLKQ IEDAGLRVAG RSGDDQLVEI          IEVNPHPWFV ACQFHPEFTS TPRDGHPLFA GFVKAASEFQ KRQAK</p>
<b>Source</b>	Baculovirus
<b>Target Names</b>	pyrG
<b>Protein Names</b>	Recommended name: CTP synthase EC= 6.3.4.2 Alternative name(s): CTP synthetase UTP--ammonia ligase
<b>Expression Region</b>	1-545
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.