



Recombinant Escherichia coli O17:K52:H18 Hydroxyacylglutathione hydrolase (gloB)

Product Code	CSB-YP010121EOG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7N874
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MNLNSIPAFD DNYIWVLNDE AGRCLIVDPG DAEPVLNAIA ANNWQPEAIF LTHHHHDHVG GVKELVKKFP QIVVYGPQET QDKGTTRVVK DGEIAFVLGH EFSVIATPGH TLGHICYFSK PYLFCGDTLF SGGCGRLFEG TPSQMYQSLK KLSALPDDTL VCCAHEYTLS NMKFALSILP HDLSINDYYR KVKELRAKNQ ITLPVILKNE RQINVFLRTE DIDLINVINE ETLQPEER FAWLRSKKDR F
Source	Yeast
Target Names	gloB
Protein Names	Recommended name: Hydroxyacylglutathione hydrolase EC= 3.1.2.6 Alternative name(s): Glyoxalase II Short name= Glx II
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.