



Recombinant Escherichia coli O45:K1 Choline dehydrogenase (betA)

Product Code	CSB-MP005335EOK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7MCD0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O45:K1 (strain S88 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	<p>MQFDYIIIIGA GSAGNVLATR LTEDPNTTVL LLEAGGPDYR FDFRTQMPAA LAFPLQGKRY NWAYETEPEP FMNRRMECG RGKGLGGSSL INGMCIYIRGN ALDLNWAQE PGLNWSYLD CLPYRKAET RDVGENDYHG GDGPVSVTTS KPGVNPLFEA MIEAGVQAGY PRTDDLNGYQ QEGFGPMDRT VTPHGRRAST ARGYLDQAKS RPNLTIRTHA MTDHIFFDGK RAVGVEWLEG DSTIPTRAAA NKEVLLCAGA IASPQILQRS GVGNAELLA EFDIPLVHELP GVGENLQDHL EMYLQYECKE PVS LYPALQW WNQPKIGAEW LFGGTGVGAS NHFEAGGFIR SREEFAWPNI QYHFLPVAIN YNGSNAVKEH GFQCHVGSMR SPSRGHVRIK SRDPHQHPAI LFNYSHEQD WQEFRDAIRI TREIMHQPAL DQYRGREISP GVECQTDEQL DEFVRNHAET AFHPCGTCKM GYDEMAVVDA EGRVHGLEGL RVVDASIMPQ IITGNLNATT IMIGEKMADM IRGKDALPRS TARYFVANGM PVRACKMSRD VN</p>
Source	Mammalian cell
Target Names	betA
Protein Names	Recommended name: Choline dehydrogenase Short name= CDH Short name= CHD EC= 1.1.99.1
Expression Region	1-562
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.