



# Recombinant Escherichia coli O7:K1 Argininosuccinate lyase (argH)

<b>Product Code</b>	CSB-EP002213EON
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B7NU40
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MALWGGRFTQ AADQRFKQFN DSLRFDYRLA EQDIVGSVAW SKALVTVGVL TAEFAQLEE ALNVLLEDVR ARPQQILESD AEDIHSWVEG KLIDKVGQLG KKLHTGRSRN DQVATDLKLW CKDTVSELLT ANRQLQSALV ETAQNNQDAV MPGYTHLQRA QPVTFAHWCL AYVEMLARDE SRLQDALKRL DVSPLGCGAL AGTAYEIDRE QLAGWLG FAS ATRNSLDSVS DRDHVLELLS AAAIGMVHLS RFAEDLIFFN TGEAGFVELS DRVTS GSSLM PQKKNPDALE LIRGKCGRVQ GALTGMMMTL KGLPLAYNKD MQEDKEGLFD ALDTWLDCLH MAALVLDGIQ VKRPRCQEAQ QQGYANATEL ADYLVAKGVP FREAHHIVGE AVVEAIRQ GK PLEDLPLSEL QKFSQVIGED VYPILSLQSC LDKRAAKGGV SPQQVAQAIA FAQARLG
<b>Source</b>	E.coli
<b>Target Names</b>	argH
<b>Protein Names</b>	Recommended name: Argininosuccinate lyase Short name= ASAL EC= 4.3.2.1 Alternative name(s): Arginosuccinase
<b>Expression Region</b>	1-457
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.