



Recombinant Danio rerio Proline-rich nuclear receptor coactivator 2 (pnrc2)

Product Code	CSB-MP018290DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7TWP7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MGGGERYNIP DRPAPKKSQP VSRGKQRSRD QNGVMHSAAS GALGVPHHLR RGEKGTYSYSW SPEARQAVSV DKKNGSVRFA TPYDQNWESH LNKLLSAQCG QNYAGAKFSE PPSPSVLPKP PSHWVSLPMG DHRELMTFQL KSLKLVQA
Source	Mammalian cell
Target Names	pnrc2
Protein Names	Recommended name: Proline-rich nuclear receptor coactivator 2
Expression Region	1-148
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.