



Recombinant Human Stress-70 protein, mitochondrial (HSPA9)

Product Code	CSB-EP010830HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P38646
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p> ASEA IKGAVVGIDL GTTNSCVAVM EGKQAKVLEN AEGARTTPSV VAFTADGERL VGMPAKRQAV TNPNTFYAT KRLIGRRYDD PEVQKDIKNV PFKIVRASNG DAWVEAHGKL YSPSQIGAFV LMKMKETAEN YLGHTAKNAV ITVPAYFNDS QRQATKDAGQ ISGLNVL RVI NEPTAAALAY GLDKSEDKVI AVYDLGGGTF DISILEIQKG VFEVKSTNGD TFLGGEDFDQ ALLRHIVKEF KRETGVDLTK DNMALQVRVRE AA EKAKCELS SSVQTDINLP YLTMDSSGPK HLNMKLTRAQ FEGIVTDLIR RTIAPCQKAM QDAEVSKSDI GEVILVGGMT RMPKVQQTVQ DLFGRAPSKA VNPDEAVAIG AAIQGGVLAG DVTDVLLLDV TPLSLGIETL GGVFTKLINR NTTIPTKKSQ VFSTAADGQT QVEIKVCQGE REMAGDNKLL GQFTLIGIPP APRGVPQIEV TFDIDANGIV HVSAKDKGTG REQQIVIQSS GGLSKDDIEN MVKNAEKYAE EDRRKKERVE AVNMAEGIIH DTETKMEEFK DQLPADECNK LKEEISKMRE LLARKDSETG ENIRQAASSL QQASLKL FEM AYKKMASERE GSGSSGTGEQ KEDQKEEKQ </p>
Source	E.coli
Target Names	HSPA9
Protein Names	Recommended name: Stress-70 protein, mitochondrial Alternative name(s): 75 kDa glucose-regulated protein Short name= GRP-75 Heat shock 70 kDa protein 9 Mortalin Short name= MOT Peptide-binding protein 74 Short name=
Expression Region	47-679
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	The product encoded by this gene belongs to the heat shock protein 70 family which contains both heat-inducible and constitutively expressed members. The latter are called heat-shock cognate proteins. This gene encodes a heat-shock cognate protein. This protein plays a role in the control of cell proliferation. It may also act as a chaperone.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.