



Recombinant Bovine Pregnancy-associated glycoprotein 1

Product Code	CSB-EP644917BO-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q29432
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	RGSNLTT HPLRNIKDLV YMGNITIGTP PQEFQVVFDT ASSDLWVPSD FCTSPACSTH VRFRHLQSST FRLTNKTFRI TYGSGRMKGV VVHDTVIRIGN LVSTDQPFGL SIEEYGFEGR IYDGVLGLNY PNISFSGAIP IFDKLKNQRA ISEPVFAFYL SKDEREGSVV MFGGVDHRYE EGELNWWPLI QAGDWSVHMD RISIERKIIA CSDGCKALVD TGTSDIVGPR RLVNNIHRLI GAIPRGSEHY VPCSEVNTLP SIVFTINGIN YPVPGRAYIL KDDRGRCYTT FQENRVSSST ETWYLGDVFL RLYFSVFDRG NDRIGLARAV
Source	E.coli
Protein Names	Recommended name: Pregnancy-associated glycoprotein 1 Short name= PAG 1 EC= 3.4.23.-Alternative name(s): Pregnancy-specific protein B Short name= PSP-B
Expression Region	54-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.