



# Recombinant Arabidopsis thaliana Cyclin-B2-2 (CYCB2-2)

<b>Product Code</b>	CSB-BP654394DOA
<b>Abbreviation</b>	CYCB2-2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q39070
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MVNPEENNRRN LVVVKPITEIL QDDDKRSRKF GVEMKRQNRRL ALGVINHNLV GAKAYPCVVN KRRGLSQRKQ ESCDKKKLDS LHPSISRSQE ETKKPKPSGN EFGDCIFIDE EEEKNEEVTL DQPMPMSLEE PYIEFDPMEE EVEMEDMEEE QEEPVLIDIDE YDANNSLAAV EYVQDLYDFY RKTERFSCVP LDYMAQQFDI SDKMRAILID WLIEVHDKFE LMNETLFLTV NLIDRFLSKQ AVARKKLQLV GLVALLLACK YEEVSVPIVE DLVVISDKAY TRTDVLEMEK IMLSTLQFNM SLPTQYPFLK RFLKAAQSDK KLEILASFLI ELALVDYEMV RYPPSLLAAT AVYTAQCTIH GFSEWNSTCE FHCHYSENQL LECCRRMVRL HQKAGTDKLT GVHRKYSSSK FGYIATKYEA AHFLVSDSH
<b>Source</b>	Baculovirus
<b>Target Names</b>	CYCB2-2
<b>Protein Names</b>	Recommended name: Cyclin-B2-2 Alternative name(s): Cyc2b-At Cyclin-2b G2/mitotic-specific cyclin-B2-2 Short name= CycB2;2
<b>Expression Region</b>	1-429
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.