



Recombinant Human 60S ribosomal protein L13a (RPL13A)

Product Code	CSB-EP020130HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P40429
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AEVQVLVLD GRGHLLGRLA AIVAKQVLLG RKVVVVRCEG INISGNFYRN KLKYLAFRLK RMNTNPSRGP YHFRAPSRIF WRTVRGMLPH KTKRGQAALD RLKVFDGIPP PYDKKKRMVV PAALKVVRLK PTRKFAYLGR LAHEVGWKYQ AVTATLEEKR KEKAKIHYRK KKQLMRLRKQ AEKNVEKKID KYTEVLKTHG LLV
Source	E.coli
Target Names	RPL13A
Protein Names	Recommended name: 60S ribosomal protein L13a Alternative name(s): 23 kDa highly basic protein
Expression Region	2-203
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L13P family of ribosomal proteins. It is located in the cytoplasm. Transcript variants utilizing alternative polyA signals have been observed. This gene is co-transcribed with the small nucleolar RNA genes U32, U33, U34, and U35, which are located in its second, fourth, fifth, and sixth introns, respectively. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.