






Recombinant Epstein-Barr virus Replication and transcription activator (BRLF1), partial

Product Code	CSB-EP662061EFC
Abbreviation	Recombinant Epstein-Barr virus BRLF1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3KSS7
Product Type	Recombinant Protein
Immunogen Species	Epstein-Barr virus (strain GD1) (HHV-4) (Human herpesvirus 4)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	PRLRAPKSRRTSRPNRGQTPCPSNAEEPEQPWIAAVHQESDERPIFPHPSKP TFLPPVKRKKGLRDSREGMFLPKPEAGSAISDVFEQREVCQPKRIRPFHPPGS PWANRPLPASLAPTPTGPVHEPVGSLTPAPVPKPLDPAPAVTPEASHLLEDPD EETSQAVKALREMA DTVIPQKEEAAICGQMDLNHPPPRGHLDELTTTLESMT DLNLDSP LTP ELNEILD TFLNDECLLHAMHISTGLSIFDTSLF
Research Area	Others
Source	E.coli
Target Names	BRLF1
Protein Names	Recommended name: Replication and transcription activator Short name= Rta Alternative name(s): Immediate-early protein Rta
Expression Region	352-605aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	34.9 kDa
Protein Length	Partial
Image	



We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.