



Recombinant Mouse ETS translocation variant 1 (Etv1)

Product Code	CSB-MP007854MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P41164
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MDGFYDQQVP YVVTNSQRGR NCTEKPTNVR KRKFINRDLA HDSEELFQDL SQLQETWLAE AQVPDNDQF VPDYQAESLA FHGLPLKIKK EPHSPCSELG SACSQEQPFK FSYGEKCLYN VSAYDQKPQV GMRPSNPPTP SSTPVSPPLHH ASPNTAHTPK PDRAFFAHL P PSQSIPDSTY PMDHRFRRQL SEPCNSFPPL PTMPREGRPM YQRQMSEPNI PFPPQGFKQE YHDPVYEHTT MVGGAASQSF PPPLMIKQEP RDFAYDSEVP SCHSIYMRQE GFLAHPSTRTE GCMFEKGPRQ FYDDTCVVPE KFDGDIKQEP GMYREGPTYQ RRGSLQLWQF LVALLDDPSN SHFIAWTGRG MEFKLEPEE VARRWGIQKN RPAMNYDKLS RSLRYYYEKG IMQKVAGERY VYKFVCDPEA LFSMAFPDNQ RPLLKTDMER HINEEDTVPL SHFDESMTYM PEGGCCNPHP YNEGYVY
Source	Mammalian cell
Target Names	Etv1
Protein Names	Recommended name: ETS translocation variant 1 Alternative name(s): Ets-related protein 81
Expression Region	1-477
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.