




Recombinant Human Keratin, type II cuticular Hb5 (KRT85)

Product Code	CSB-EP012580HU
Abbreviation	Recombinant Human KRT85 protein
Uniprot No.	P78386
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.201.
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MSCRSYRISSGCGVTRNFSSCSAVAPKTGNRCCISAAPYRGVSCYRGLTGFG SRSLCNLGSCGPRIAVGGFRAGSCGRSFGYRSGGVCGPSPPCITTTSVNESL LTPLNLEIDPNAQCQVKQEEKEQIKSLNSRFAAFIDKVRFLQEQNKLLQTKWQFY QNQRCCESNLEPLFSGYIETLRREAECVEADSGRLASELNHVQEVLEGYKKKY EEEVALRATAENEFVVLKDVDCAYLRKSDLEANVEALVEESSFLRRLYEEEIR VLQAHISDTSVIVKMDNSRDLNMDCIAEIKAQYDDVASRSRAEAEESWYRSKCE EMKATVIRHGETLRRTKEEINELNRMQRLTAEIENAKCQRAKLEAAVAEAEQQ GEAALSDARCKLAELEGALQKAKQDMACLLKEYQEVMSKLGLEDIEIATYRRL EGEEHRLCEGVGSVNVCSVSSSRGGVSCGGLSYSTTPGRQITSGPSAIGGSITV VAPDSCAPCQPRSSSFSCGSSRSVRFA
Research Area	Cell Biology
Source	E.coli
Target Names	KRT85
Expression Region	1-507aa
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	62.7 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.