



Recombinant Eubacterium rectale Uracil-DNA glycosylase (ung)

Product Code	CSB-MP025641EOX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C4Z8Z2
Product Type	Recombinant Protein
Immunogen Species	Agathobacter rectalis (strain ATCC 33656 / DSM 3377 / JCM 17463 / KCTC 5835 / VPI 0990) (Eubacterium rectale)
Purity	>85% (SDS-PAGE)
Sequence	MSMLDGEWAK AVGAEFKKPY YRDLYNFVRD EYATHVVYPP ADDIFNAMHF TPLDKVKVLI LGQDPYHNVN QAHGLSFSVL PSQKDIPPSL QNIYKELHDD LGCDIPNNGY LKKWADQGV LLLNTVLTVRA HQANSHQGHG WEKFTDAIIE AVNEQDRPIV YMLWGRPAQS KIPMLTNPKE LILKAPHPSP LSAYRGFFGC KHFSKANEFL KEHGVEPVDW QIEDI
Source	Mammalian cell
Target Names	ung
Protein Names	Recommended name: Uracil-DNA glycosylase Short name= UDG EC= 3.2.2.27
Expression Region	1-225
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.