



Recombinant Human Cytochrome b-c1 complex subunit 2, mitochondrial (UQCRC2)

Product Code	CSB-MP025668HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P22695
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>SLKVAP KVKATAAPAG APPQPQDLEF TKLPNGLVIA SLENYSPVSR IGLFIKAGSR YEDFSNLGTT HLLRLTSSLT TKGASSFKIT RGIEAVGGKL SVTATRENMA YTVECLRGDV DILMEFLNV TTAPEFRRWE VADLQPQLKI DKAVAFQNPQ THVIENLHAA AYRNALANPL YCPDYRIGKV TSEELHYFVQ NHFTSARMAL IGLGVSHPV LKQVAEQFLNM RGGLGLSGAK ANYRGGEIRE QNGDSLVAHA FVAESAVAGS AEANAFSVLQ HVLGAGPHVK RGSNTTSHLH QAVAKATQQP FDVSAFNASY SDSGLFGIYT ISQATAAGDV IKAAYNQVKT IAQGNSNTD VQAAKNKLKA GYLMSVESSE CFLEEVGSQA LVAGSYMPPS TVLQQIDSV NADIINAACK FVSGQKSMAS SGNLGHTPFV DEL</p>
Source	Mammalian cell
Target Names	UQCRC2
Protein Names	<p>Recommended name: Cytochrome b-c1 complex subunit 2, mitochondrial Alternative name(s): Complex III subunit 2 Core protein II Ubiquinol-cytochrome-c reductase complex core protein 2</p>
Expression Region	15-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>