



Recombinant Bovine NADH dehydrogenase [ubiquinone] iron-sulfur protein 3, mitochondrial (NDUFS3)

Product Code	CSB-MP015662BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P23709
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	ES SAADTRPTVR PRNDVAHKQL SAFGEYVAEI LPKYVQQVQV SCFNELEICI HPDGVIPVLT FLRDHSNAQF KSLADLTAVD IPTRQNRFEI VYNLLSLRFN SRIRVKTYTD ELTPIESSVP VYKAANWYER EIWDMFGVFF ANHPDLRRIL TDYGFEGHPF RKDFPLSGYV ELRYDDEVKR VVAEPVELAQ EFRKFDLNSP WEAFPAYRQP PESLKLEAGD TKPEAK
Source	Mammalian cell
Target Names	NDUFS3
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] iron-sulfur protein 3, mitochondrial EC= 1.6.5.3 EC= 1.6.99.3 Alternative name(s): Complex I-30kD Short name= CI-30kD NADH-ubiquinone oxidoreductase 30 kDa subunit
Expression Region	39-266
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.