




Recombinant Rickettsia conorii Outer membrane protein A (ompA), partial

Product Code	CSB-EP706579RMS
Relevance	Elicits protective immunity.
Abbreviation	Recombinant Rickettsia conorii ompA protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q52657
Product Type	Recombinant Protein
Immunogen Species	Rickettsia conorii (strain ATCC VR-613 / Malish 7)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	DMDAKFGAWISPFVGNATQKMCNSISGYKSDTTGGTIGFDGFVSDDLVLGLAY TRADTDIKLKNNKTGDKNKVESNIYSLYGLYSVPYENLFVEAIASYS DNKIRSKS RRVIATTLETVG YQTANGKYKSESYTGQLMAGYTYMMSENINLTPLAGLRYSTI KDKSYKETGTTYQNLTVKGKNYNTFDGLLGAKVSSNINVNEIVLTPELYAMVDY AFKNKVSAIDARLQGMTAPLPTNSFKQSKTSFDVGVGVTA KHKMMMEYGINYDT NIGSKYFAQQGSVKVRVNF
Research Area	Microbiology
Source	E.coli
Target Names	ompA
Protein Names	190 kDa antigen Cell surface antigen rOmp A
Expression Region	1734-2021aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged
Mol. Weight	51.8kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.