



Recombinant *Cricetulus griseus* Cyclin-dependent kinase 2 (CDK2)

Product Code	CSB-BP005061DXU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O55076
Product Type	Recombinant Protein
Immunogen Species	<i>Cricetulus griseus</i> (Chinese hamster) (<i>Cricetulus barabensis griseus</i>)
Purity	>85% (SDS-PAGE)
Sequence	MENFQKVEKI GEGTYGVVYK AKNKLTGEVV ALKKIRLDTE TEGVPSTAIR EISLLKELNH PNIVKLLDVI HTENKLYLVF EFLHQDLKKF MDASAVTGIP LPLIKSYLFQ LLQGLAFCHS HRVLHRDLKP QNLLINAEGS IKLADFGLAR AFGVPVRTYT HEVVTWLWYRA PEILLGCKYY STAVDIWSLG CIFAEMVTRR ALFPGDSEID QLFRIFRITLG TPDEVVWPGV TSMPDYKPSF PKWARQDFSK VVPPLDEDGR SLLSQMLHYD PNKRISAKAA LAHPFFQDVT KPVPHLRL
Source	Baculovirus
Target Names	CDK2
Protein Names	Recommended name: Cyclin-dependent kinase 2 EC= 2.7.11.22 Alternative name(s): Cell division protein kinase 2
Expression Region	1-298
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.