



Recombinant Human Bone morphogenetic protein 15 (BMP15)

Product Code	CSB-BP002735HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95972
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	QAD GISAEVTASS SKHSGPENNQ CSLHPFQISF RQLGWDHWII APPFYTPNYC KGTCLRVLDR GLNSPNHAI QNLINQLVDQ SVPRPSCVPY KYVPISVLM I EANGSILYKE YEGMIAESCT CR
Source	Baculovirus
Target Names	BMP15
Protein Names	Recommended name: Bone morphogenetic protein 15 Short name= BMP-15 Alternative name(s): Growth/differentiation factor 9B Short name= GDF-9B
Expression Region	268-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.