



Recombinant Arabidopsis thaliana Adenylate kinase 1 (ADK1)

Product Code	CSB-MP001508DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O82514
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MATGGAAADL EDVQTVDLMS ELLRRLKCSQ KPDKRLIFIG PPGSGKGTQS PVVKDEYCLC HLSTGDMLRA AVASKTPLGV KAKEAMEKGE LVSDDLVVGI IDEAMNKPKC QKGFILDGFP RTVTQAEKLD EMLKRRGTEI DKVLNFAIDD AILEERITGR WIHPSSGRSY HTKFAPPKTP GVDDITGEPL IQRKDDNADV LKSRLAAFHS QTQPVIDYYA KKAULTNIQA EKAPQEVTSE VKKALS
Source	Mammalian cell
Target Names	ADK1
Protein Names	Recommended name: Adenylate kinase 1 Short name= AK 1 EC= 2.7.4.3 Alternative name(s): ATP-AMP transphosphorylase 1
Expression Region	1-246
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.