



# Recombinant Mouse Lymphocyte antigen 6E (Ly6e)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP717533MO-B   |
| <b>Abbreviation</b>      | Ly6e   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q64253   |
| <b>Storage Buffer</b>    | Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0   |
| <b>Product Type</b>      | Recombinant Proteins   |
| <b>Immunogen Species</b> | Mus musculus (Mouse)   |
| <b>Purity</b>            | Greater than 85% as determined by SDS-PAGE.  |
| <b>Sequence</b>          | LMCFSCTDQK NNINCLWPVS CQEKDHYCIT LSAAAGFGNV NLGYTLNKGC<br>SPICPSENVN LNLGVASVNS YCCQSSFCNF SA  |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | Name:Ly6eSynonyms:Ly67, Sca-2, Tsa-1   |
| <b>Protein Names</b>     | Recommended name: Lymphocyte antigen 6E Short name= Ly-6EAlternative name(s): Stem cell antigen 2 Thymic shared antigen 1 Short name= TSA-1  |
| <b>Expression Region</b> | 27-108   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |