




# Recombinant Mouse All-trans-retinol dehydrogenase [NAD (+)] ADH7 (Adh7)

<b>Product Code</b>	CSB-EP717544MO
<b>Relevance</b>	Could function in retinol oxidation for the synthesis of retinoic acid, a hormone important for cellular differentiation.
<b>Abbreviation</b>	Recombinant Mouse Adh7 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q64437
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MGTAGKVIKCKAAVLWGVNQPFSEIEVAPPKAKEVRVKILATGICRTDDHIIK GSMVSKFPVIVGHEAVGVVESVGEVTTVRPGDKVIPLFLPQCRCNACLNPE GNLCIRSDLTGRGVLADGTTTRFTCKGKPVQHFMNTSTFTEYTVLDESSVAKVD GAAPPEKACLIGCGFSTGYGAAVKTAKVTPGSTCVVFLGGVGLSVIMGCKAA GASRIIGIDINKDKFQKALAVGATECISPKDSTKPISEVLSDMTGNTIQYTFEVI RLETMVDALSSCHMNYGTSVVVGAPPSAKMLTYDPMLLFTGRTWKGCVFVGG WKSRRDDVPKLVTEFLEKKFDLDQLITHLTPFNNINEGFELLYSGKSIRTVLTF
<b>Research Area</b>	Metabolism
<b>Source</b>	E.coli
<b>Target Names</b>	Adh7
<b>Expression Region</b>	1-374aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	55.9kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.