



Recombinant Mouse NAD (P)H dehydrogenase [quinone] 1 (Nqo1)

Product Code	CSB-MP717562MO
Abbreviation	Nqo1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64669
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	AARRALIVL AHSEKTSFNY AMKEAAVEAL KKRGWEVLES DLYAMNFNPI ISRNDITGEL KDSKNFQYPS ESSLAYKEGR LSPDIVAEHK KLEAADLVIF QFPLQWFGVP AILKGWFERV LVAGFAYTYA AMYDNGPFQN KKTLLSITG GSGSMYSLQG VHGD MNVILW PIQSGILRFC GFQVLEPQLV YSIGHTPPDA RMQILEGWKK RLETVWEETP LYFAPSSLFD LNFQAGFLMK KEVQEEQKKN KFGLSVGHHL GKSIPADNQI KARK
Source	Mammalian cell
Target Names	Nqo1
Protein Names	Recommended name: NAD(P)H dehydrogenase [quinone] 1 EC= 1.6.5.2 Alternative name(s): Azoreductase DT-diaphorase Short name= DTD Menadione reductase NAD(P)H:quinone oxidoreductase 1 Phylloquinone reductase Quinone
Expression Region	2-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene is a member of the NAD(P)H dehydrogenase (quinone) family and encodes a cytoplasmic 2-electron reductase. This FAD-binding protein forms homodimers and reduces quinones to hydroquinones. This protein's enzymatic activity prevents the one electron reduction of quinones that results in the production of radical species. Mutations in this gene have been associated with tardive dyskinesia (TD), an increased risk of hematotoxicity after exposure to benzene, and susceptibility to various forms of cancer. Altered expression of this protein has been seen in many tumors and is also associated with Alzheimer's disease (AD). Alternate transcriptional splice variants, encoding different isoforms, have been characterized.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.