



Recombinant Human Fc receptor-like protein 6 (FCRL6), partial

Product Code	CSB-BP718779HU
Abbreviation	Recombinant Human FCRL6 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DN72
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	LYLQAWPNPVFEGDALTLRCQGWNKNTPLSQVKFYRDGKFLHFSKENQTLISM GAATVQSRGQYSCSGQVMYIPQTFTQTSETAMVQVQELFPPPVLSAIPSPEPR EGSLVTLRCQTKLHPLRSALRLLFSFHKDGHTLQDRGPHPELCIPGAKEGDSG LYWCEVAPEGGQVQKQSPQLEVRVQAPVSRPVLTLHHGPADPAVGDMVQLL CEAQRGSPPIYSFYLDEKIVGNHSAPCGGTTSLLPVKSEQDAGNYSCEAEN SVSRERSEPKKLSLKGSQVLFTPASNW
Research Area	Immunology
Source	Baculovirus
Target Names	FCRL6
Protein Names	Recommended name: Fc receptor-like protein 6 Short name= FcR-like protein 6 Short name= FCRL6 Alternative name(s): Fc receptor homolog 6 Short name= FcRH6 IFGP6
Expression Region	20-307aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	34.2 kDa
Protein Length	Extracellular Domain
Image	

