



# Recombinant Human Moesin (MSN)

<b>Product Code</b>	CSB-EP015048HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P26038
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>PKTISVRVT TMDAELEFAI QPNTTGKQLF DQVVKTIGLR EVWFFGLQYQ  DTKGFSTWLK LNKKVTAQDV RKESPLLKFK RAKFY PEDVS EELIQDITQR  LFFLQVKEGI LNDDIYCPPE TAVLLASYAV QSKYGDFNKE VHKSGYLAGD  KLLPQRVLEQ HKLNKDQWEE RIQVWHEEHR GMLREDAVLE YLKIAQDLEM  YGVNYFSIKN KKGSELWLGV DALGLNIYEQ NDRLTPKIGF PWSEIRNISF  NDKKFVIKPI DKKAPDFVfy APRLRINKRI LALCMGNHEL YMRRRKPDTI  EVQQMKAQAR EEKHQKQMER AMLENEKKKR EMAEKEKEKI  EREKEELMER LKQIEEQTKK AQQELEEQTR RALELEQERK RAQSEAEKLA  KERQEAEAAK EALLQASRDQ KKTQEQLALE MAELTARISQ LEMARQKKES  EAVEWQQKAQ MVQEDLEKTR AELKTAMSTP HVAEPAENEQ  DEQDENGAEA SADLRADAMA KDRSEEERTT EAEKNERVQK HLKALTSELA  NARDESKKTA NDMIHAENMR LGRDKYKTLR QIRQGNTKQR IDEFESM</p>
<b>Source</b>	E.coli
<b>Target Names</b>	MSN
<b>Protein Names</b>	Recommended name: Moesin Alternative name(s): Membrane-organizing extension spike protein
<b>Expression Region</b>	2-577
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	Moesin (for membrane-organizing extension spike protein) is a member of the ERM family which includes ezrin and radixin. ERM proteins appear to function as cross-linkers between plasma membranes and actin-based cytoskeletons. Moesin is localized to filopodia and other membranous protrusions that are important for cell-cell recognition and signaling and for cell movement.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.