



Recombinant Escherichia coli Uroporphyrinogen decarboxylase (hemE)

Product Code	CSB-EP025678ENP-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B1IUQ0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Purity	>85% (SDS-PAGE)
Sequence	MTELKNDRYL RALLRQPVDV TPVWMMRQAG RYLPEYKATR AQAGDFMSLC KNAELACEVT LQPLRRYPLD AAILFSDILT VPDAMGLGLY FEAGEGPRFT SPVTCKADVD KLPIPDPEDE LGYVMNAVRT IRRELKGEVP LIGFSGSPWT LATYMVEGGS SKAFTVIKKM MYADPQALHA LLDKLAHSV LYLNAQIKAG AQAVMIFDTW GGVLTGRDYQ QFSLYYMHKI VDGLLRENDG RRVPVTLFTK GGGQWLEAMA ETGCDALGLD WTTDIADARR RVGNKVALQG NMDPSMLYAP PARIEEEVAT ILAGFGHGEG HVFNLGHGHIH QDVPPEHAGV FVEAVHRLSE QYHR
Source	E.coli
Target Names	hemE
Protein Names	Recommended name: Uroporphyrinogen decarboxylase Short name= UPD Short name= URO-D EC= 4.1.1.37
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.