



Recombinant Escherichia coli 6-phosphogluconolactonase (pgl)

Product Code	CSB-BP017861ENP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B1IXL9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Purity	≥85% (SDS-PAGE)
Sequence	MKQTVYIASP ESQQIHVWNL NHEGALTLTQ VVDVPGQVQP MVVSPDKRYL YVGV RPEFRV LAYRIAPDDG ALTFAAESAL PGSPHISTD HQGQFVFGS YNAGNVSVTR LEDGLPVGWV DVVEGLDGCH SANISPDNRT LWVPALKQDR ICLFTVSDDG HLVAQDPAEV TTVGAGPRH MVFHPNEQYA YCVNELNSSV DVWELKDPHG NIECVQTLDM MPENFSDTRW AADIHITPDG RHLIACDRTA SLITVFSVSE DGSVLSKEGF QPTETQPRGF NVDHSGKYLI AAGQKSHHIS VYEIVGEQGL LHEKGRYAVG QGPMWVVVNA H
Source	Baculovirus
Target Names	pgl
Protein Names	Recommended name: 6-phosphogluconolactonase Short name= 6-P-gluconolactonase EC= 3.1.1.31
Expression Region	1-331
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.