



Recombinant Mouse Guanine nucleotide-binding protein subunit alpha-13 (Gna13)

Product Code	CSB-MP009585MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P27601
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MADFLPSRSV LSVCFPGCVL TNGEAEQQRK SKEIDKCLSR EKTYVKRLVK ILLGAGESG KSTFLKQMRI IHGQDFDQRA REEFRPTIYS NVIKGMRLV DAREKLHIPW GDNKNQLHGD KLMAFDTRAP MAAQGMVETR VFLQYLPAIR ALWEDSGIQN AYDRRREFQL GESVKYFLDN LDKLGVPDYI PSQQDILLAR RPTKGIHEYD FEIKNVPFKM VDVGGQRSER KRWFECDSDV TSILFLVSSS EFDQVLMEDR QTNRLTESLN IFETIVNNRV FSNVSIILFL NKTDLLEEKV QVVSIDYFL EFEGDPHCLR DVQKFLVECF RGKRRDQQQR PLYHHFTTAI NTENIRLVFR DVKDTILHDN LKQLMLQ
Source	Mammalian cell
Target Names	Gna13
Protein Names	Recommended name: Guanine nucleotide-binding protein subunit alpha-13 Short name= G alpha-13 Short name= G-protein subunit alpha-13
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.