



# Recombinant Bovine F-actin-capping protein subunit beta (CAPZB)

<b>Product Code</b>	CSB-EP004514BO-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P79136
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MHPCRRSLPF PLNCQLVKVG TADYGGASDQ SDQQLDCALD LMRRLPPQQI EKNLSIDLIDL VPSLCEDLLS SVDQPLKIAR DKVVGKDYLL CDYNRDGDSY RSPWSNKYDP PLEDGAMPSA RLRKLEVEAN NAFDQYRDLY FEGGVSSVYL WDLDHGFAGV ILIKKAGDGS KKIKGWDSI HVVEVQEKSS GRTAHYKLT TVMLWLQTNK SGSGTMNLGG SLTRQMEKDE TVSDCSPHIA NIGRLVEDME NKIRSTLNEI YFGKTKDIVN GLRSVQTFAD KSKQEALKND LVEALKRKQQ C
<b>Source</b>	E.coli
<b>Target Names</b>	CAPZB
<b>Protein Names</b>	Recommended name: F-actin-capping protein subunit beta Alternative name(s): CapZ beta
<b>Expression Region</b>	1-301
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	CAPZB is a member of the F-actin capping protein family. This gene encodes the beta subunit of the barbed-end actin binding protein. The protein regulates growth of the actin filament by capping the barbed end of growing actin filaments.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.