



Recombinant *Aeromonas hydrophila* subsp. *hydrophila* Uracil-DNA glycosylase (ung)

Product Code	CSB-BP025641AUH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A0KMZ1
Product Type	Recombinant Protein
Immunogen Species	<i>Aeromonas hydrophila</i> subsp. <i>hydrophila</i> (strain ATCC 7966 / DSM 30187 / JCM 1027 / KCTC 2358 / NCIMB 9240)
Purity	>85% (SDS-PAGE)
Sequence	MQTWTDVIGS EKEQDYFKAT LATVRSEREA GKVIYPPATD VFNAFKLTEL DQVKVVILGQ DPYHGPNQAH GLCFSVLPGV RTPPSLVNIY KEMQRDLPGF VTPSHGFLES WAKQGVLLLN TVLTVQAGQA HSHAHLGWET FTDRVIEQIN ASCQGVVFLW GAHAQKKGR FIDRQRHHVL SAPHSPLSA HRGFIGCGHF SETNRLLSQQ GMSPINWHSV CE
Source	Baculovirus
Target Names	ung
Protein Names	Recommended name: Uracil-DNA glycosylase Short name= UDG EC= 3.2.2.27
Expression Region	1-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.