



ADORA2A Antibody, FITC conjugated

Product Code	CSB-PA001376LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29274
Immunogen	Recombinant Human Adenosine receptor A2a protein (291-412AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adenosine receptor A2a, ADORA2A, ADORA2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ADORA2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.