



AICDA Antibody, HRP conjugated

Product Code	CSB-PA001487LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZX7
Immunogen	Recombinant Human Single-stranded DNA cytosine deaminase protein (1-60AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Single-stranded DNA cytosine deaminase (EC 3.5.4.38) (Activation-induced cytidine deaminase) (AID) (Cytidine aminohydrolase), AICDA, AID
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	AICDA