





ARHGAP10 Antibody

Product Code	CSB-PA002012LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A1A4S6
Immunogen	Recombinant Human Rho GTPase-activating protein 10 protein (542-744AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho GTPase-activating protein 10 (GTPase regulator associated with focal adhesion kinase 2) (Graf-related protein 2) (Rho-type GTPase-activating protein 10), ARHGAP10, GRAF2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ARHGAP10

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA002012LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA002012LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.