



AURKA Antibody, FITC conjugated

Product Code	CSB-PA002454HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14965
Immunogen	Recombinant Human Aurora kinase A protein (1-131AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aurora kinase A (EC 2.7.11.1) (Aurora 2) (Aurora/IPL1-related kinase 1) (ARK-1) (Aurora-related kinase 1) (hARK1) (Breast tumor-amplified kinase) (Serine/threonine-protein kinase 15) (Serine/threonine-protein kinase 6) (Serine/threonine-protein kinase aurora-A), AURKA, AIK AIRK1 ARK1 AURA AYK1 BTAK IAK1 STK15 STK6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	AURKA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.