




BUB1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA002882LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O43683 |
| Immunogen | Recombinant Human Mitotic checkpoint serine/threonine-protein kinase BUB1 protein (256-577AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, IHC; Recommended dilution: IHC:1:20-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Mitotic checkpoint serine/threonine-protein kinase BUB1 (hBUB1) (EC 2.7.11.1) (BUB1A), BUB1, BUB1L |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | BUB1 |

Image



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA002882LA01HU at dilution of 1:100