




BUB1 Antibody

Product Code	CSB-PA002882LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43683
Immunogen	Recombinant Human Mitotic checkpoint serine/threonine-protein kinase BUB1 protein (256-577AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitotic checkpoint serine/threonine-protein kinase BUB1 (hBUB1) (EC 2.7.11.1) (BUB1A), BUB1, BUB1L
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BUB1

Image



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA002882LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.