



BUB1B Antibody, HRP conjugated

Product Code	CSB-PA002883LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60566
Immunogen	Recombinant Human Mitotic checkpoint serine/threonine-protein kinase BUB1 beta protein (230-394AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitotic checkpoint serine/threonine-protein kinase BUB1 beta (EC 2.7.11.1) (MAD3/BUB1-related protein kinase) (hBUBR1) (Mitotic checkpoint kinase MAD3L) (Protein SSK1), BUB1B, BUBR1 MAD3L SSK1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BUB1B