



C9 Antibody, Biotin conjugated

Product Code	CSB-PA004259LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02748
Immunogen	Recombinant Human Complement component C9 protein (74-155AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement component C9 [Cleaved into: Complement component C9a; Complement component C9b], C9
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	C9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.