



CD83 Antibody, FITC conjugated

Product Code	CSB-PA004962LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01151
Immunogen	Recombinant Human CD83 antigen protein (20-144AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD83 antigen (hCD83) (B-cell activation protein) (Cell surface protein HB15) (CD antigen CD83), CD83
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD83
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.