




CLCN1 Antibody

Product Code	CSB-PA005480LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35523
Immunogen	Recombinant Human Chloride channel protein 1 protein (634-988AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chloride channel protein 1 (CIC-1) (Chloride channel protein, skeletal muscle), CLCN1, CLC1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CLCN1

Image



Immunofluorescent analysis of HepG2 cells using CSB-PA005480LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)