





CNOT1 Antibody

Product Code	CSB-PA005663LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A5YKK6
Immunogen	Recombinant Human CCR4-NOT transcription complex subunit 1 protein (604-867AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CCR4-NOT transcription complex subunit 1 (CCR4-associated factor 1) (Negative regulator of transcription subunit 1 homolog) (NOT1H) (hNOT1), CNOT1, CDC39 KIAA1007 NOT1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CNOT1

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA005663LA01HU at dilution of 1:100



Immunofluorescent analysis of U251 cells using CSB-PA005663LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.