



CXCR3 Antibody, FITC conjugated

Product Code	CSB-PA006253LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49682
Immunogen	Recombinant Human C-X-C chemokine receptor type 3 protein (4-50AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-X-C chemokine receptor type 3 (CXC-R3) (CXCR-3) (CKR-L2) (G protein-coupled receptor 9) (Interferon-inducible protein 10 receptor) (IP-10 receptor) (CD antigen CD183), CXCR3, GPR9
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CXCR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.