



DNAJA1 Antibody, Biotin conjugated

Product Code	CSB-PA006999LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31689
Immunogen	Recombinant Human DnaJ homolog subfamily A member 1 protein (8-106AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DnaJ homolog subfamily A member 1 (DnaJ protein homolog 2) (HSDJ) (Heat shock 40 kDa protein 4) (Heat shock protein J2) (HSJ-2) (Human DnaJ protein 2) (hDj-2), DNAJA1, DNAJ2 HDJ2 HSJ2 HSPF4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DNAJA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.