



DNMT1 Antibody, Biotin conjugated

Product Code	CSB-PA007082HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26358
Immunogen	Recombinant Human DNA (cytosine-5)-methyltransferase 1 protein (1288-1632AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA (cytosine-5)-methyltransferase 1 (Dnmt1) (EC 2.1.1.37) (CXXC-type zinc finger protein 9) (DNA methyltransferase Hsal) (DNA MTase Hsal) (M.Hsal) (MCMT), DNMT1, AIM CXXC9 DNMT
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DNMT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.