



EPM2A Antibody, Biotin conjugated

Product Code	CSB-PA007738LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B3EWF7
Immunogen	Recombinant Human Laforin, isoform 9 protein (1-331AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Laforin, isoform 9, EPM2A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	EPM2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.