



F5 Antibody, Biotin conjugated

Product Code	CSB-PA007929HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12259
Immunogen	Recombinant Human Coagulation factor V protein (1490-1614AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coagulation factor V (Activated protein C cofactor) (Proaccelerin, labile factor) [Cleaved into: Coagulation factor V heavy chain; Coagulation factor V light chain], F5
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	F5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.