





FABP7 Antibody

Product Code	CSB-PA007956LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15540
Immunogen	Recombinant Human Fatty acid-binding protein, brain protein (2-132AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fatty acid-binding protein, brain (Brain lipid-binding protein) (BLBP) (Brain-type fatty acid-binding protein) (B-FABP) (Fatty acid-binding protein 7) (Mammary-derived growth inhibitor related), FABP7, BLBP FABPB MRG
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	FABP7


Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA007956LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human glioma using CSB-PA007956LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA007956LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.