



FABP7 Antibody, Biotin conjugated

Product Code	CSB-PA007956LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15540
Immunogen	Recombinant Human Fatty acid-binding protein, brain protein (2-132AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fatty acid-binding protein, brain (Brain lipid-binding protein) (BLBP) (Brain-type fatty acid-binding protein) (B-FABP) (Fatty acid-binding protein 7) (Mammary-derived growth inhibitor related), FABP7, BLBP FABPB MRG
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	FABP7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.