





RIPOR2 Antibody

Product Code	CSB-PA008285LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4F9
Immunogen	Recombinant Human Rho family-interacting cell polarization regulator 2 protein (1-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho family-interacting cell polarization regulator 2, RIPOR2, C6orf32 DIFF48 FAM65B KIAA0386 PL48
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RIPOR2

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA008285LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA008285LA01HU at dilution of 1:100