



RIPOR2 Antibody, Biotin conjugated

Product Code	CSB-PA008285LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4F9
Immunogen	Recombinant Human Rho family-interacting cell polarization regulator 2 protein (1-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho family-interacting cell polarization regulator 2, RIPOR2, C6orf32 DIFF48 FAM65B KIAA0386 PL48
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RIPOR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.